



Protect eyes and face against liquids and dust with LATITUDE™ Protective Face Shield. Manufactured in the USA with 3M material, our full face shield offers maximum protection on all sides of the face, protecting the eyes, nose and mouth from contact with liquids and dust. LATITUDE™ is not impact resistant nor intended for use against open flames or sparks.

Together we are addressing not only the needs of the medical community on the front-line fighting COVID-19 but also the needs of industry to support their reopening during the current pandemic.



**LIGHTWEIGHT & OPTICALLY CLEAR**

Provides clear, unobstructed view. Lightweight and clear. Effective for use against liquids and dust. 3M Shield Technology includes Anti-Fog coating.



**MAXIMUM PROTECTION**

Complete coverage including the sides of the face preventing all contact of liquids and dust from eyes, mouth, and nose.



**COMFORTABLE AND FUNCTIONAL**

LATITUDE™ comes with a comfortable black foam. It also easily goes over glasses or goggles. Latex, Fiberglass and Distortion Free.

**LATITUDE - PROTECTIVE FACE SHIELD**

**Description**

**LATITUDE™ PLUS GENERIC HEADBAND**  
3M™ Anti-Fog Face Shield Film, 9960FS, 6.8 mil

**LATITUDE™ PLUS CUSTOM HEADBAND**  
3M™ Anti-Fog Face Shield Film, 9960FS, 6.8 mil  
Custom 4-color image with logo

**Anti Fog    Per Bag    Per Case    MOQ**



1

200

200



1


200

4800

**HEADBAND OPTIONS**

**LATITUDE™ PLUS:**

Black Label - Black Foam/Elastic



**LATITUDE™ PLUS Custom:**

4 Color Custom Artwork



**LATITUDE™ PLUS - CUSTOM HEADBANDS**

Increase brand awareness while keeping your patrons safe.

Design Guidelines: <http://design.graphiccontrols.com>

**1 Shield / Bag  
200 Shields / Case**



8 5/8"

12 3/8"

3-13-01-074 (4/2020)